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Instruction for use of DURASPLINT[®] LC

The transparent light-curing resin
specifically designed to adjust
DURAN[®] splints.

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We recommend the above mentioned items for the fabrication of adjusted occlusal splints made of DURAN® to facilitate your work and to ensure a user-friendly production process.

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Instructions regarding work safety

Avoid skin contact with uncured resins. Make sure to wear gloves during processing. Wear protective clothing (smock, protective goggles, safety mask). Switch on suction unit.

Instructions regarding light polymerisation

Light polymerisation requires a wavelength between 350 nm and 450 nm. The LC-6 Light Oven offers this spectral range.

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Use flat trimmed models. Mount models in an articulator following the dentist's instructions and using the recommended bite registration. Optionally, you can use magnetic systems.

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Measure the model using the parallelometer taking into account the direction of insertion and mark the prosthetic equator.

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Mark splint extension in the area of slight undercuts in red.

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We recommend blocking out strong undercuts, like space closures, using BLUE-BLOKKER® or SIL-KITT.

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The same applies for interdental spaces and bridge elements.

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Pressurise an ISO-FOLAN® foil first before thermoforming DURAN® 1.5 mm or 2.0 mm according to the instructions.

We recommend using the pellet cover.

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Important: It is essential to use DURAN® as basis for the light-curing material.

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Reduce the pressure-moulded DURAN® roughly by means of the HM carbide bur or the diamond separating disc.

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Hint: For optimum fit remove splint only after light-polymerisation of DURASPLINT® LC and after reducing the splint.

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The material surface in the complete splint extension should be sandblasted using aluminium oxide (Al_2O_3) or roughened with a DIMO® wheel (brown, coarse).

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Take off grease residues from roughened area with ethyl alcohol to guarantee a good compound. Apply one thin layer of LC-Primer on the complete splint extension and proceed with curing for 5 minutes in LC-6 Light Oven.

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Apply LC Insulation to the antagonist model with a brush, to prevent DURASPLINT® LC material from sticking to the model.

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Adapt one bar of the DURASPLINT® material on the pressure moulded splint. Additionally, CLEAR-BLOKKER® can be applied as modelling aid to reach a smooth result.

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Important: Avoid skin contact with uncured material. Wear gloves while processing!

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For easier modelling moisten your fingertips with LC Insulation. Now you can model the material to the desired shape.

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Hint: We recommend applying a thin layer of the material, slightly exceeding the splint extension, in order to avoid visible material transitions between DURASPLINT® LC and DURAN®.

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Insert movement patterns and impressions using the articulator. Movements should be made only once in the corresponding direction before returning to the centric position with the articulator being open.

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Hint: Repeated movements or return movements with the articulator being closed might cause distortion due to the slightly sticky material.

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Remove the model with the pressure moulded splint from articulator and place it centrally in the LC-6 Light Oven. Then cure the material 2 times for 10 minutes. If the LC-6 Light Oven is used for a serial workflow, please make sure to cool down the LC-6 Light Oven after each polymerisation period by opening the drawer.

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After light-curing you can finish the splint still being on the working model.

The slightly sticky inhibition layer can be removed by using isopropyl alcohol.

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For initial finishing we recommend using LC Burs A, B or C

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After grinding in sliding surfaces and teeth cusps according to the therapeutic instructions you can proceed with finishing the splint.

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Reduce material up to the splint extension carefully by means of the HM Carbide Cutter or the diamond separating disc. Avoid any contact between carbide cutter and model.

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Remove splint carefully from the model, taking it off simultaneously from both sides.

Hint: This is essential to avoid cracks that might occur if the splint is removed from one side only.

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A chamfered POLYFIX disc is best suited for interdental finishing of the splint.

For further finishing we recommend using the LC burs A, B and C.

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
We recommend using the DIMO® disc coarse to smooth the surfaces of the rims or optionally sandpaper (both items are not included in the DURASPLINT® LC kit).

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The SD Polishing Brush, SD-Felt Wheel and SD-High-Gloss Polishing Mop are best suited for polishing in combination with SD-Polishing Paste A and B.

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Polished DURAN® splint with DURASPLINT® LC build-up on DURAN®.

Hint: We do not recommend occlusal polishing, in order to maintain the full functionality of the adjusted splint surfaces.

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You can use CLEAR-BLOKKER® in a layer of 1-2 mm for any kind of completions or corrections. Make sure the inhibition layer stays on the object. If this layer has already been removed, apply LC-Primer to the corresponding area and cure for 5 minutes. Then apply CLEAR-BLOKKER® directly from the syringe using an application cannula. Proceed with light-curing for 10 minutes, e. g. in the LC-6 Light Oven.

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Additional information:

DURASPLINT® LC and CLEAR-BLOKKER® do not contain methyl methacrylate monomer (only the LC-Primer contains methyl methacrylate monomer). Please use only light curing ovens that do not develop high temperature levels during operation (we recommend the LC-6 Light Oven).

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The CETRON® Care Set is the best choice for daily cleansing and care.

CE

DE/GB 500/11/18 G REF PM0094.01