IMPRIMO® LC Model

Produktinformation

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IMPRIMO® LC Model

Product information

Description:	IMPRIMO® LC Model is a light-curing, methacrylate-based resin that is specifically engineered for the fabrication of dental models. Photo-polymerisation is done by means of image projections systems (385 nm). The printed models are characterized by similar haptic and optic qualities as conventional plaster models. IMPRIMO® LC Model should be used in building processes with a z-resolution of 10-100 µm. Printing parameters available for Asiga MAX™ and Asiga PRO 4K.
Processing instructions:	Make sure IMPRIMO® LC Model is tempered to 23C – 30° C prior to use. Fill the material into the tray wearing your personal protection equipment (safety glasses, gloves and face mask). Small air bubbles can be removed carefully with a spatula. In case of contamination of the material such as dirt or fragments, the material can be filtrated and repeatedly mixed due to its low viscosity. Leave the material to stand for approx. one hour to avoid any air inclusions.
Contra-indications:	IMPRIMO® LC Model should not be used for any purpose other than the generative CAD/CAM fabrication of dental models and is not suitable for oral application. IMPRIMO® LC Model is a laboratory product.
Post-processing:	Upon completion of the building process, immediate further processing is recommended. If this cannot be guaranteed, the produced objects should be left in the liquid IMPRIMO® Model resin. Unpolymerised resin residues can be removed completely using Isopropanol and a soft brush or may be resolved in an appropriate cleaning device such as IMPRIMO® Clean or RS wash using IMPRIMO® Cleaning Liquid.
Post-curing:	Final and complete curing of the objects is done in an appropriate curing unit such as IMPRIMO® Cure or RS cure. The use of protective gas or of a vacuum to prevent the creation of a dispersion layer is crucial. A wavelength of 385 nm is required. Polymerisation in the IMPRIMO® Cure unit takes 5 minutes with a with a 20-minute nitrogen flush. In the RS cure unit polymerization takes 5 minutes in a vacuum environment. The corresponding parameters are stored in the internal programmes of the RS cure and the IMPRIMO® Cure units. Any deviation from the described post-curing process can lead to modifications of the mechanical properties or the material colour.
Shelf life / storage:	Store the material at room temperature (18 °C to 28 °C) and protect it from light. Any exposure to light might damage the material. Carefully seal the container immediately after use.
Hazard warnings:	Prolonged skin contact with un-cured material and inhalation of monomer vapours can lead to irritations. In single cases, allergic reactions to certain ingredients of IMPRIMO® LC Model might occur, for more details, please refer to our safety data sheets on IMPRIMO® LC Model. In the event of eye contact, rinse thoroughly with water and seek medical advice. In the event of skin contact, rinse thoroughly with running water and soap. Attention: Polymerised resins are chemically resistant. Avoid stains on clothing
Disposal:	According to local authority prescriptions. Cured objects can be disposed of with other household waste. Unpolymerized material has to be disposed of at a waste collection point. All information on processing our material - verbal, written or practical - is given to the best of our knowledge and must be read as references without obligation. Material use and processing is beyond our control and lies within the responsibility of the user. For professional use only.
Technical data:	Density: approx. 1.1 g/m Viscosity: approx. 0.7 Pa s Cured material (depends on irradiation unit): Flexural strength: 108 MPa Elongation at break: 5 % Shore hardness: 85 D <small>*These data result from measurements of a representative sample.</small> Contains: alkoxyliated bisphenol-A dimethacrylate, initiators, stabilizers, dyes and pigments. Material Safety Data Sheet: www.scheu-dental.com/en/service

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